



CERTIFICATE OF APPROVAL OF A POLYESTER RESIN

Certificate No. **MATS/2437/4**

This certificate is issued to the company named below. The resin described has been examined in accordance with the requirements of Lloyd's Register and is approved for use in constructions built under Lloyd's Register's survey. This approval is subject to Lloyd's Register being informed of any changes in or modifications to the resin and the product being used in accordance with the manufacturer's instructions and with the relevant requirements of Lloyd's Register's Rules and Regulations.

Company	CARLO RICCO & F.LLI SPA CORREGGIO ITALY
Trade name	POLIPLAST
Resin	LR 96
Application	Laminating resin
Type	Dicyclopentadiene
Characteristics	
Applicable LR Rules	Rules and Regulations for Classification of Special Service Craft
Approved Variants	See attached Appendix

Valid until **1 August 2021**

Date **8 February 2017**

M Jogia
Lead Specialist to Lloyd's Register EMEA
A member of the Lloyd's Register group



CERTIFICATE OF APPROVAL A OF POLYESTER RESIN

Appendix to certificate MATS/2437/4

Approved Variants

LR 96 X15	Pre-accelerated, Thixotropic, Short geltime
LR 96 X25	Pre-accelerated, Thixotropic, Medium geltime
LR 96 X45	Pre-accelerated, Thixotropic, Long geltime
LR 96 X W	Pre-accelerated, Thixotropic, Medium geltime, low styrene emission
LR 96 X KB15	Pre-accelerated,,Thixotropic, Short geltime, low exotherm peak
LR 96 X KB25	Pre-accelerated,,Thixotropic, Medium gel time, low exotherm peak
LR 96 X KB45	Pre-accelerated,,Thixotropic, Long geltime, low exotherm peak
LR 96 X KB W	Pre-accelerated,,Thixotropic, Low exotherm peak, low styrene emission
LR 96 XD15	Pre-accelerated,,Thixotropic, Short geltime, low viscosity
LR 96 XD25	Pre-accelerated,,Thixotropic, Medium geltime, low viscosity
LR 96 XD45	Pre-accelerated,,Thixotropic, Long geltime , low viscosity
LR 96 XD W	Pre-accelerated,,Thixotropic, Short - medium - long geltime, low styrene emission
LR 96 XD KB15	Pre-accelerated,,Thixotropic, Short geltime, low exotherm peak
LR 96 XD KB25	Pre-accelerated,,Thixotropic, Medium geltime, low exotherm peak
LR 96 XD KB45	Pre-accelerated,,Thixotropic, Long geltime, low exotherm peak
LR 96 XD KB W	Pre-accelerated,,Thixotropic, Short - medium - long geltime, low exotherm peak, low styrene emission
LR 96.02	Not pre-accelerated, low viscosity, for infusion
LR 96.02 D	Not pre-accelerated, low viscosity, for infusion
LR 96 INF	Not pre-accelerated, low viscosity, for infusion
LR 96.06 INF	Pre-accelerated, Low viscosity, for infusion

Date **8 February 2017**

M Jorgia

Lead Specialist to Lloyd's Register EMEA
A member of the Lloyd's Register group